LANAFORM® Filo® LCD Screen Thermometer

within 10 minutes.
- The «Filo» pacifier thermometer is designed to work for 200 hours. The «Filo» symbol will flash if the battery is low.
- When this happens, simply replace your pacifier with a new «Filo» pacifier thermometer, ensuring you comply with the regulations concerning waste disposal.
- If for any reason you need to interrupt the measurement, press the procedure described above.

NOTES:
- Oral temperature measurements are usually 0.5 °C lower than the actual body temperature.
- If the body temperature is above 42.9 °C, the LCD screen will display «Hi».
MAINTENANCE AND CLEANING
- You can sterilize the thermometer by immersing the silicon latex teat ONLY in warm water.
Please note:
- You must NEVER completely immerse the pacifier thermometer in boiling water, as you risk damaging the electronic components.
- Do not use abrasive detergents or chemical solvents.
ADVICE REGARDING THE DISPOSAL OF WASTE
- All packaging is composed of materials that pose no hazard for the environment and which can be disposed of in your local sorting centre to be used as secondary raw material. The cardboard box will be disposed of in a paper recycling bin. The packaging film must be taken to your local sorting and recycling centre.
- When you have finished using the device, please dispose of it in an environmentally friendly way and in accordance with the laws. Prior to disposal, please remove the battery and dispose of it at a collection point that it can be recycled.

EC 0029 DECLARATION OF CONFORMITY
Medical device - Category Ia.
We, the company LANAFORM, declare that the medical device «Filo» LA090110 complies with Council Directive 93/42/EEC concerning medical devices and its most recent amendments 2007/45/EC, APFRAGZ, 156 Chaussée de Willevoorde, 1120 Brussels (Belgium), notified body no. 0029, has checked this product (according to Annex IV).

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PLEASE READ ALL THE INSTRUCTIONS BEFORE USING YOUR «PACIFIER THERMOMETER», ESPECIALLY THIS BASIC SAFETY INFORMATION:
- Do not place the thermometer in direct sunlight, or expose it to high temperatures, humidity or dirt.
- Do not use this device in a room in which aerosols (sprays) are used or in a room in which oxygen is administered.
- Do not immerse the device in water.
- Do not drop the device and do not use it if it is damaged.
- This device must be used in an ambient temperature between 10°C and 40°C.
- Do not expose the thermometer to electric shocks.
- Do not expose this thermometer to temperatures of greater than 60°C or below -20°C.
- Do not use the device if the relative humidity is above 90% or below 20%.
- Do not use the device in the vicinity of large electromagnetic fields, such as microwaves or portable telephones.
- It is recommended to place the thermometer in the room for around 30 minutes before use.
- Use of this thermometer is not a substitute for consulting your doctor.
- If there is a problem with your thermometer, please contact your reseller.
- Do not try to repair the thermometer yourself.
- This unit is not intended for use by people (including children) with reduced physical, sensory or mental capabilities, or persons with lack of experience or knowledge, unless a person responsible for their safety can supervise or instruct them beforehand on how to use the device. Children should be responsible to make sure they do not play with the device.
- Only use the thermometer for its intended purpose, as described in the manual.
- Before each use, examine the pacifier, checking for any cracks, tears, or other signs of use. Do not use the pacifier thermometer if there is the slightest sign of damage.
- Do not try to pull the pacifier or teat to take it out of the device or pull off the rim of the button, as you might damage the temperature sensor and electronic components, and consequently the sensitivity of the thermometer.
- Do not let children play with the device. It is not a toy. The thermometer is a sensitive electronic device which could be damaged if it is not used correctly in the presence of an adult.
- Do not attach a ribbon or string to the pacifier thermometer as your child could be strangled by it.
- Only use this thermometer for temperature measurements taken orally.
- Do not use the pacifier thermometer like a normal thermometer.
- The battery is not replaceable. Do not try to remove it.

FEATURES
- LCD screen to display the temperature.
- Last measurement in memory.
- Sound alert to announce that the maximum temperature has been reached.
- Double sound alert when the temperature exceeds 37.8°C.
- Automatic stop function after 10 minutes.
- Protective cap to transport the product hygienically
DESCRIPTION OF THE THERMOMETER:
1 - On/Off button
2 - LCD screen
3 - Pacifier
4 - Sensor
5 - Protective cap
INSTRUCTIONS FOR USE
The following table is a guide showing average human body temperatures according to age.
Normal temperatures according to age.

Normal temperatures according to age	Temp. °C	Normal conditions of use	Measurement range	Precision	Automatic stop
0 to 2 years	36.4 to 38.0	Display resolution 0.1°C	32°C to 42.9°C	±0.1°C 35-38°C ±0.2°C 38-42.0°C	10 mins
3 to 10 years	36.1 to 37.8	Operating temperature 10 to 40°C			
11 to 65 years	35.9 to 37.6	Storage temperature 20°C to 60°C			
>65 years	35.8 to 37.5	Humidity level 20%-90% Weight +/- 14 grams			

M.B.: The temperature of the human body varies throughout the day. It can also be influenced by a number of external factors: age, gender, type and thickness of skin etc.

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